

## Mounting Flange

The Mounting Flange allows the chamber for the CRC®-25R to be mounted in a cabinet or fume hood flush with the surface of the counter top, saving counter space. This is a MUST for all hot labs with limited counter space.

This mounting unit is designed to be used with the chamber for the CRC®-25R. It can also be used for the CRC®-15 Ultra, and the CRC®-127R. The flange is made of attractive brushed aluminum that lies flat against the counter top. It allows the well liner insert to be installed from the top and permits removal of the well for cleaning or decontamination.

The flush Mounting Flange is easily installed in existing or new cabinets or hoods. Drawings are available upon request.



**Mounting Flange** ..... 7310-2307