

Anthropomorphic 3D Skull Phantom





Model 603

For Rapid Assessment of Image Displacement in Gamma Knife and Other Treatment Planning Systems

Features:

- Images well on T1, T2 and 3DTOF MRI acquisitions
- Images well on CT scans
- Stereotactic frame can be applied to special reinforced pads (included)
- Images can be imported into stereotactic localization program
- CT scans can be used to assess MRI accuracy



Three dimensional orthogonal acrylic rod matrix through cranial volume enables assessment of image distortions