

Mammography Artifact Evaluation Phantom



Models 014C & 014E

Quickly detect artifacts

The American College of Radiology and MQSA recommend a uniform 4 cm thick "high grade" cassette sized phantom for evaluation of mammography artifacts as it is often difficult to identify artifacts based on clinical or standard phantom images.

CIRS has designed two phantom sets to meet these recommendations.

Each set contains 2 slabs. The Model 014C slabs each measure 18 X 24 X 2 cm thick. The Model 014E slabs each measure 24 X 30 X 2 cm thick.