

TO ABC

RADIOGRAPHY & FLUOROSCOPY

A routine test object designed to be used quickly and easily on a regular basis (e.g. weekly/monthly) to provide an ongoing check of beam alignment and centring.

Used for radiography or fluoroscopy, enables the following checks to be made:

The correspondence of the light field with the effective radiation field.


The determination of the position of the central beam.

The central beam angle deviation from 90° , relative to the plane of the film/detector, is below a limit that beyond which the image quality would be adversely affected.




For more information or to find your regional distributor please visit our website

www.leadstestobjects.com | tel. +44(0)1423 321102



TO ABC

RADIOGRAPHY & FLUOROSCOPY



For more information or to find your regional distributor please visit our website

www.leadstestobjects.com | tel. +44(0)1423 321102

Copyright © 2011 Leeds Test Objects Ltd | SpecSh 169b 07/07/2011



TO ABC

RADIOGRAPHY & FLUOROSCOPY



A special variant is available for Siemens' Vertex system. Please request at time of order.

VERTICAL BUCKY JIG

For more information or to find your regional distributor please visit our website

www.leadstestobjects.com | tel. +44(0)1423 321102