


BT-CS-137 Brachytherapy Source Holder

Positions Cs-137 Seeds in Dose Calibrator

Specially designed to fit Capintec dose calibrators, the BT-Cs-137 brachytherapy source holder positions Cs-137 seeds in the ionization chamber for accurate assays of activity in any Capintec Dose Calibrator.



Centimeter markings allow convenient positioning and determination of source length. The centrally positioned tube is easily adjusted with a thumb screw.

Cs-137 Seed Holder, Standard Length Chamber 0976-0004