

## BT-IR-192 Brachytherapy Ribbon Source Holder

Positions Ir-192 Seeds in Dose Calibrator

Specially designed to fit Capintec standard dose calibrators, the BT-Ir-192 brachytherapy ribbon-source holder positions Ir-192 seeds in the ionization chamber for accurate assays of activity.

The coil holds up to 12 Ir-192 seeds with 1-cm spacing and no overlap.



**Ir-192 Ribbon Holder, Standard Length Chamber..... 0976-0005**

**Ir-192 Ribbon Holder, BT Chamber..... 0976-0007**