

# Prostate Demonstration Phantom



**Model 058**

*Non-disposable urethane phantom  
for prostate imaging*

The Model 058 is a derivative of the Model 053 Ultrasound Prostate Training Phantom. The Model 058 is durable and appropriate for repetitive demonstration scanning.

The prostate along with structures simulating the rectal wall, seminal vesicles and urethra is contained within an 11.5 cm X 7.0 cm X 9.5 cm clear container.

A needle is embedded within the prostate to demonstrate needle localization. The phantom also contains a simulated lesion and calcification cluster.

This phantom is an ideal demonstration device for rectal scanning.

