

TGF Ruler



Model 033

Create permanent alignment record

Confirmation of alignment between laser lights and the X-ray energy field is a routine QA check. Although there is equipment available, they do not provide a permanent record of the results. The Model 033 TGF Ruler is a low cost, simple to use tool that can provide this permanent record.

The acrylic ruler has radiopaque indices at 1 mm increments making it suitable for mammography applications. The cut-out enables the TGF ruler to be compatible with the Unfors DXR+ or alone.

To use, position the zero index at the edge of the light field, expose the film and then visually inspect the visible measurement index on the film and record your results.