



# Beam Alignment Test Instrument

## GAMMEX 430

The Gammex 430 is an instrument used to analyze the alignment of a linear accelerator. It is recommended that beam alignment tests be completed at least once a year to determine the following problem situations:

- Displaced focal spot
- Asymmetric collimators
- Non-intersection of the collimator and gantry axes
- Lack of gantry arm support

The Gammex 430 complies with the alignment requirements of AAPM Task Group 40.

## SPECIFICATIONS

Construction . . . . . Acrylic with lead blocks in a  
10x10 cm field  
Size . . . . . 13x23.5x30 cm (5x9.3x11.5 in)  
Fits . . . . . 8x10 in standard film pack  
Weight . . . . . 6 kg (13.5 lbs)

