



## TOR DEN Plus

DENTAL

Leeds Test Objects' TOR DEN Plus phantom is designed to be used to check the image quality performance of modern high-resolution digital dental intra-oral x-ray systems. These tests should be carried out at acceptance to establish a baseline, and on a routine basis, keeping an on-going record of the test results, which will indicate any deterioration in performance.

TOR DEN digital is used for units with digital image acquisition. It goes beyond the requirements of standards IEC 61223-3-4 and IEC 61223-2-7 and enables the following tests to be made:

- Limiting Spatial Resolution  
(8.0 - 24 LP/mm)
- Low Contrast Resolution  
(0.2, 0.5, 1.0 and 2.0mm dia.)  
(2.0, 1.0, 0.5, 0.3 and 0.1 mm depth)
- Image Receptor Dose

TOR DEN Plus is supplied in a tough case with a bespoke foam insert.




For more information or to find your regional distributor please visit our website



# TOR DEN Plus

DENTAL



Radiograph of TOR DEN Plus

For more information or to find your regional distributor please visit our website

[www.leadstestobjects.com](http://www.leadstestobjects.com) | tel. +44(0)1423 321102